



February 2004

Engineering US1 Braid Equivalence Through Construction Scaling

Executive Summary

This development programme focused on reproducing the dimensional and structural characteristics of an established US1 braided construction using alternative fibre architectures and controlled braid scaling.

Multiple UHMWPE and polyester braid configurations, together with supplementary nylon and silk variants, were evaluated to determine how construction variables influenced braid diameter, density, coverage, and mechanical behaviour.

By systematically optimising carrier configuration, yarn selection, braid density, and process parameters, the programme established repeatable manufacturing routes capable of achieving US1 dimensional equivalence across multiple architectures.

The Challenge

Maintaining an established USP braid size while introducing alternative fibre configurations requires precise control of braid architecture.

Even small changes in carrier count, yarn denier or braid density influence:

- Finished diameter
- Surface coverage
- Braid compaction
- Mechanical stability
- Handling characteristics
- Manufacturing repeatability

The challenge was to reproduce the dimensional profile of an existing US1 braid while allowing flexibility in fibre selection and construction design.

Engineering Objectives

The programme aimed to:

- Achieving US1 dimensional equivalence across multiple braid constructions
- Optimise braid architecture through controlled construction scaling
- Evaluate the influence of carrier configuration and fibre selection
- Maintain structural stability and mechanical performance
- Establish repeatable manufacturing parameters suitable for production
- Increase flexibility for future material selection and product development

Development Approach

Construction Scaling

A series of braid architectures was developed to investigate how structural variables affected dimensional performance.

Development variables included:

- Carrier count
- Yarn linear density
- Fibre distribution
- Cover-to-core balance
- Pick count (PPI)
- Haul-off conditions
- Twist orientation

Alternative constructions combined UHMWPE with polyester, alongside selected nylon and silk variants, to evaluate

the effect of fibre composition on braid geometry and manufacturability.

Engineering Evaluation

Multiple braid configurations were produced and compared against the established US1 reference construction.

Engineering assessment focused on understanding how controlled changes in braid architecture could maintain equivalent diameter while preserving structural integrity and manufacturing stability.

Validation & Performance Assessment

Each construction was evaluated against the dimensional and mechanical requirements associated with a US1 braid.

Assessment included:

- Diameter verification
- Tensile performance
- Extension under load
- Pick count (PPI)
- Braid uniformity
- Manufacturing stability

Comparative evaluation identified braid architectures capable of achieving dimensional equivalence while maintaining consistent manufacturing performance.

Process Control & Repeatability

Consistent manufacturing performance was achieved through careful control of key braiding parameters.

Critical process controls included:

- Balanced carrier tension
- Standardised gear configurations
- Controlled haul-off conditions
- Defined twist orientation
- Uniform yarn distribution
- Continuous monitoring of braid density (PPI)

The controls minimised process variability and enabled reliable comparison across multiple braid architectures.

Engineering Outcome

The programme demonstrated that US1 braid characteristics could be successfully replicated through controlled construction scaling rather than reliance on a single fibre configuration.

Key engineering outcomes included:

- Successful replication of US1 braid dimensions
- Optimised carrier architecture for multiple fibre combinations
- Improved understanding of the relationship between fibre selection and braid geometry
- Repeatable manufacturing parameters for alternative constructions
- Increased flexibility for future material substitution and product customisation

The project reinforced the principle that braid architecture is the primary driver of dimensional performance, enabling equivalent structures to be engineered using different material combinations.

Application Relevance

The engineering principles demonstrated are directly applicable to medical braided structures where dimensional precision, mechanical consistency, and manufacturing repeatability are critical.

Surgical Sutures & Closure Systems

Braided sutures requiring precise USP size control, consistent handling, and reliable manufacturing performance.

Orthopaedics & Sports Medicine

High-strength braided fixation systems where controlled diameter and structural stability support device performance.

Cardiovascular Devices

Braided reinforcement structures and catheter components requiring precise dimensional control and repeatable braid architecture.

Minimally Invasive Surgical Devices

Compact braided components used in delivery systems and surgical instrumentation, where braid geometry directly influences performance.

Advanced Implantable Textile Systems

Hybrid braided textile platforms requiring controlled architecture, material flexibility, and scalable manufacturing processes.

Key Takeaways

- Braid architecture can be scaled to reproduce established USP size profiles across alternative fibre constructions.
- Dimensional equivalence depends on the interaction of carrier configuration, fibre distribution, and braid density rather than material selection alone.
- Construction scaling provides greater flexibility for future material substitution and product development.
- Careful process control enables consistent manufacture across multiple braid architectures.
- Precision braid engineering supports repeatable, production-ready medical textile development.

“Controlled construction scaling enabled multiple fibre architectures to achieve the dimensional characteristics of an established US1 braid without compromising manufacturing repeatability.”



Contact Corza Medical today

