



July 2006

Tapered Braided Constructs

Executive Summary

A medical device OEM required the development of a tapered braided construction transitioning from US2 to US0 for applications requiring controlled diameter reduction without sacrificing strength, structural integrity, or handling performance.

Rather than treating as a simple geometric change, the programme focused on engineering a braid architecture capable of maintaining consistent functional behaviour across changing dimensions.

Through optimisation of braid construction, material selection, and process control, multiple tapered configurations were successfully developed and validated, demonstrating repeatable production of complex variable-diameter braided structures.

The Challenge

Producing a tapered braid introduces multiple interacting variables that become increasingly difficult to control as diameter decreases.

The smallest diameter region (US0) represented the highest risk area due to increased sensitivity to structural instability and localised performance loss.

Key engineering challenges included:

- Maintaining braid integrity during controlled fibre removal
- Retaining mechanical strength at reduced diameters
- Preventing abrupt stiffness transitions across taper regions
- Managing different tension and deformation responses between UHMWPE and polyester.
- Achieving precise, repeatable taper geometry
- Ensuring scalability and manufacturing consistency

Successfully resolving these competing requirements required simultaneously control of structure, materials, and process.

Engineering Objectives

The programme focused on:

- Developing tapered braid constructions transitioning from US2 to US0
- Maintaining structural continuity throughout taper regions
- Preserving tensile performance and handling characteristics
- Evaluating hybrid and polyester-only construction strategies
- Establishing repeatable manufacturing conditions
- Creating a scalable framework for variable-diameter braid production

Development Approach

Three braid architectures were developed to compare the influence of material selection and taper geometry on performance and manufacturability:

- Hybrid UHMWPE/polyester single taper
- Hybrid UHMWPE/polyester double taper
- Polyester single taper

Each design used controlled fibre removal and defined braid architecture to create smooth stiffness transitions while maintaining structural integrity and repeatable manufacture.

Process Optimisation

Repeatable taper formation depended on precise control of braid geometry throughout production.

Optimisation focused on:

- Controlled fibre removal and reintegration
- Carrier balance and yarn tension management
- Clock-number sequencing to define taper geometry
- Monitoring braid density through transition zones
- Controlled braid angle and finished diameter

Validated process settings were established for each construction, enabling taper geometry to become a controlled engineering parameter rather than a source of manufacturing variation.

Validation & Performance Assessment

Constructs were evaluated to confirm both mechanical performance and geometric consistency.

Assessment included:

- Break load evaluation, particularly within the US0 region
- Runnage consistency
- Extension behaviour
- Diameter verification across taper sections

- Picks per inch (PPI)
- Taper uniformity and transition quality

Testing confirmed that controlled tapering could be achieved while maintaining structural continuity and repeatable mechanical performance across all configurations.

Process Control & Repeatability

Repeatable taper formation was achieved through careful control of the manufacturing process and braid architecture.

Key process controls included:

- Controlled fibre removal and reintegration sequences
- Carrier balance and yarn tension management
- Clock-number sequencing to define taper geometry
- Monitoring of braid angle, diameter, and transition quality
- Standardised process parameters for each taper configuration

These controls established a repeatable manufacturing process, enabling consistent production of variable-diameter braided constructs while minimising batch-to-batch variation.

Engineering Outcome

The programme established a repeatable engineering approach for producing variable-diameter braided constructs with controlled mechanical transitions.

Key engineering outcomes included:

- Controlled taper formation through optimised braid architecture and fibre removal
- Maintained mechanical continuity across changing diameters
- Defined performance trade-offs between hybrid and mono-material constructions
- A scalable manufacturing process for repeatable variable-diameter braid production

Application Relevance

Orthopaedics & Sports Medicine

Variable-diameter fixation constructs for ligament repair, tendon reconstruction, and soft tissue stabilisation.

Surgical Sutures & Closure Systems

Braided structures requiring controlled profile transitions while maintaining handling and tensile performance.

Advanced Textile Implants

Engineered textile systems requiring region-specific mechanical behaviour and controlled load transfer.

Key Takeaways

- Tapering is not a geometric adjustment; it is a controlled structural engineering process.
- Mechanical performance across changing diameters depends on coordinated management of braid architecture and manufacturing conditions.
- Controlled fibre removal and process synchronisation improve transition stability.
- Hybrid and single-material constructions offer different balances of performance and manufacturability.
- Repeatable taper formation enables scalable production of advanced variable-diameter braided structures.

“By controlling geometry, material selection, and process inputs simultaneously, variable-diameter braids can behave as a single continuous structure.”



Contact Corza Medical today

