

SEE THE POWER OF PRECISION.

Sharppoint MICROSUTURES

Microsuture Catalog



Operating as one. Focused on you.


Surgical Specialties Corporation, a trusted source of medical devices for over 40 years, is now part of Corza Medical, a new medical device company that unites Surgical Specialties Corporation with the TachoSil hemostatic patch from Takeda Pharmaceuticals to create a leading provider of advanced surgical products and technologies.

Corza Medical's mission is to aid and champion the surgical community with outstanding products and remarkable service that save time, money and effort. We provide medical devices as well as biosurgery products and technologies that serve surgeons and their teams, distributors and OEM partners around the world.

With an expanding platform of industry-leading brands, including Surgical Specialties precision products, our goal is to deliver an exceptional customer experience and outstanding value to the entire surgical community.



Non-Absorbable

Material	Material Description	Qualities
NYLON™  • ETHILON™ * Monosof™/ Dermalon™	Monofilament <ul style="list-style-type: none"> • Dyed (black) • Packaged dry or in wetting solution Common uses: General & Skin Closure, Plastic Surgery, Ophthalmology, Microsurgery, Cardiovascular Surgery, Gastrointestinal	<ul style="list-style-type: none"> • Passes easily through tissue

• Comparable Ethicon™ Material * Comparable Covidien™ Material


Needle Specification

Needle Code	Wire Diameter	Point Style	Curvature (Degrees)	Needle Length (mm)	Chord Length (mm)	Radius (mm)
DR4	50 Micron	Taper	135	4.0	3.1	1.7
VRM4	70 Micron	M.E.T.	105	4.0	3.5	2.2
DRM4	70 Micron	M.E.T.	135	4.0	3.1	1.7
DRM5	70 Micron	M.E.T.	135	5.0	3.9	2.1
BRM5	70 Micron	M.E.T.	95/107	4.8	2.9	.8/1.9
DR4	70 Micron	Taper	135	4.0	3.1	1.7
VRM5	100 Micron	M.E.T.	105	5.0	4.3	2.7
DRM5	100 Micron	M.E.T.	135	5.0	3.9	2.1
HRM4	100 Micron	M.E.T.	150	4.5	3.3	1.7
BRM5	100 Micron	M.E.T.	95/107	4.8	2.9	.8/1.9
VR5	100 Micron	Taper	105	5.0	4.3	2.7
DR5	100 Micron	Taper	135	5.0	3.9	2.1
HSV6	100 Micron	Vas Cutting	149	6.4	4.7	2.4
VRM6	140 Micron	M.E.T.	105	6.0	5.2	3.3
DRM6	140 Micron	M.E.T.	135	6.0	4.7	2.5
VR5	140 Micron	Taper	105	5.0	4.3	2.7
DR6	140 Micron	Taper	135	6.0	4.7	2.5

Needle Technology

M.E.T.[™] (Micro Edge Taper) Needle











Micro Edge Taper technology provides surgeons with a proprietary pyramid point with exceptional sharpness for positive tissue pick-up and effortless penetration.



- Available in sizes 8-0 through 11-0
- Variety of curvatures: ¼", ⅜", and ½" circle

Nylon

M.E.T™ Needle Point



				Size			
Needle	Wire Diameter	Material	Length	11-0	10-0	9-0	8-0
VRM4 4mm 1/4 Circle 	70 Microns	Black Nylon	4"/10cm		AA-0101		
VRM5 5mm 1/4 Circle 	100 Microns	Black Nylon	5"/13cm		AA-0106		
VRM6 6mm 1/4 Circle 	140 Microns	Black Nylon	5"/13cm			AA-0116	
DRM4 4mm 3/8 Circle 	70 Microns	Black Nylon	4"/10cm	AA-0102	AA-0103		
			5"/13cm			AA-2586	
DRM4 4mm 3/8 Circle 	70 Microns	Black Nylon	1"/2.5cm		AA-2507		
DRM5 5mm 3/8 Circle 	70 Microns	Black Nylon	5"/13cm		AA-2547		
			2.5"/6.5cm			AA-2628	
			5"/13cm		AA-0111	AA-0112	
DRM5 5mm 3/8 Circle 	100 Microns	Black Nylon	1"/2.5cm			AA-2544	
DRM6 6mm 3/8 Circle 	140 Microns	Black Nylon	5"/13cm			AA-0121	AA-0122
HRM4 4mm 1/2 Circle 	100 Microns	Black Nylon	5"/13cm			AA-0105	
BRM5 5mm Bi-curve 	70 Microns	Black Nylon	1"/2.5cm		AA-2492		
			1"/2.5cm		AA-1824		
			6"/15cm		AA-2334		

Nylon

Taper Point

				Size			
Needle	Wire Diameter	Material	Length	11-0	10-0	9-0	8-0
VR5 5mm 1/4 Circle 	100 Microns	Black Nylon	5"/13cm			AA-0137	
	140 Microns		5"/13cm				AA-0140
DR4 4mm 3/8 Circle 	50 Microns	Black Nylon	4"/10cm	AA-0130	AA-0131		
	70 Microns		4"/10cm		AA-0133		
DR5 5mm 3/8 Circle 	100 Microns	Black Nylon	5"/13cm		AA-0141	AA-0142	
DR6 6mm 3/8 Circle 	140 Microns	Black Nylon	5"/13cm				AA-0145

Vas Cutting Needles

				Size			
Needle	Wire Diameter	Material	Length	11-0	10-0	9-0	8-0
HSV6 6mm 1/2 Circle 	100 Microns	Black Nylon	5"/13cm			AA-1825	
HSV6 6mm 1/2 Circle 	100 Microns	Black Nylon	1"/2.5cm			AA-2545	

Microsurgical Sutures, Lab Packs & Training Materials

As part of our commitment to microsurgical training facilities, we offer a variety of training tools to allow surgeons to practice and develop their microsurgical techniques. The training manual and laboratory support products allow the novice to develop technical competence in a laboratory setting.

Laboratory Packs-Non-Sterile Sutures



M.E.T.™ Needle Point

Needle	Wire Diameter	Material	Length	Size			
				11-0	10-0	9-0	8-0
DRM4 4mm 3/8 Circle 	70 Microns	Black Nylon	5"/13cm			AK-0100	

Taper Point

Needle	Wire Diameter	Material	Length	Size			
				11-0	10-0	9-0	8-0
VR5 5mm 1/4 Circle 	140 Microns	Black Nylon	5"/13cm				AK-2508
DR4 4mm 3/8 Circle 	50 Microns	Black Nylon	5"/13cm	AK-2119	AK-2120		
	70 Microns				AK-0105		
DR5 5mm 3/8 Circle 	100 Microns	Black Nylon	5"/13cm		AK-0106	AK-0107	
DR6 6mm 3/8 Circle 	140 Microns	Black Nylon	5"/13cm			AK-0109	AK-2546

Laboratory Manual & PracticePak™

Product	Description	Part Number
 Manual	"A Laboratory Manual for Microvascular and Microtubal Surgery"	82-0116
 PracticePak™	Latex membrane in frame (This product contains natural rubber latex which may cause allergic reactions)	AK-9000



Order Online at www.corza.com

Order by Phone 877-991-1110

Sharppoint™ Microsuture is also available through major healthcare distributors.

For more information, visit corza.com

Corza Medical, 247 Station Drive, Suite NE1, Westwood, MA 02090, U.S.A.
Phone: 610.404.3520 • 877.991.1110 (USA and Canada only) • Fax: 610.404.4010
All trademarks are property of their respective owners. ©2021 Corza Medical.
All Rights Reserved. COR-049 R3 09/21