

**blink** SINGLE-USE  
INSTRUMENTS

# Single-use Instruments and packs



# Welcome

Blink Medical is one of the world's largest manufacturers of quality, sterile, single use surgical instruments and bespoke instrument sets and packs. We supply instruments from our UK headquarters, to both NHS and private sector hospitals for use in either theatres or out patient departments. Through our exclusive distributor network, we export to over 20 countries worldwide including Australia, New Zealand and throughout continental Europe.

Offering the highest quality instruments, we provide products for Ophthalmic, Dermatology, Plastic Surgery, Cosmetic Surgery and ENT departments. More recently we have also introduced a wider range of Bi Polar Forceps and we have increased the range of standard off the shelf packs.

Since launching in 2004, our primary objective has been to supply the highest quality surgical instruments, in conjunction with first class customer service. Our high standards and attention to detail are reflected throughout the company, beginning with our manufacturing sites which adhere to ISO13485:2016, through to quality inspection, packaging, sterilisation, and the onward supply to the end user.

If you require any further information please call our customer service team on +44 121 386 8433. Alternatively to view our extensive product range visit us at [www.blinkmedical.com](http://www.blinkmedical.com).







# blink

SINGLE-USE  
INSTRUMENTS



[www.blinkmedical.com](http://www.blinkmedical.com)

## Manufacturing Facility

All Blink Medical's surgical instruments are manufactured to the highest standard. Each of our manufacturing sites are ISO accredited and audited regularly to ensure compliance. All instruments undergo stringent inspection and cleaning protocols before being shipped to our state of the art distribution centre.



## Logistics

Our UK Logistics Centre handles thousands of shipments every year, sending consignments to the UK and around the world. Whether ordering one box of forceps or a large consignment of customised sets, our highly trained professional staff are here to ensure your shipment is expertly handled and delivered on time.



## Customer Service

Our customer service department will handle your order professionally and courteously, expediting it in a timely fashion. Should you have a query or require a quotation, our friendly, helpful staff will aim to deal with your request politely and quickly. Where possible we will process any quotation for a bespoke set the same day.



## UK Quality Control

Our UK based quality control department is staffed by fully trained and experienced associates, who undertake the inspection of all instruments in accordance with the requirements of our ISO13485:2016 accredited Quality Management System.

All instruments are subject to stringent inspection in accordance with their intended use in a safe and effective manner. Only instruments meeting the acceptance criteria can be subsequently released for final packaging operations.

All instruments and procedure packs are assembled within a validated 'Class 8' cleanroom, with subsequent sterilisation of the packaged products completed ensuring sterility is achieved and maintained in accordance with the relevant regulatory requirements.

All production processes are subject to appropriate quality assurance control procedures ensuring that all product and regulatory requirements are met on a perpetual basis.



## Internationally Available

Our exclusive International Distributor network covers the majority of Europe, Australia and New Zealand. Our main distributors carry the Ophthalmic range of instruments, however some also now carry the Cosmetic, ENT and Dermatology lines as well. Should you represent a company interested in becoming a Blink Distributor or alternatively if you are a potential customer wishing to locate your nearest Blink distributor, please visit our website at [www.blinkmedical.com](http://www.blinkmedical.com) and click on the appropriate tab.

# General Surgery

## Standard Sets

### 095 - Punch Biopsy Set



Qty	Code	Description
1	HR121	Forceps Adson Fine Toothed
1	HR125	Forceps Adson Non Toothed 1.5mm
1	HR521	Needle Holder Halsey Smooth
1	20-7002	Scissors Iris Straight 10cm Sharp Tip
5		Swab Non Woven 10x10cm
1		Tray with Gallipot Holder
1		Wrap

### 185 - Biopsy Set Full



Qty	Code	Description
1	HR519	Halsey Needle Holder Smooth Narrow Jaw
1	HR122	Forceps Adson Toothed 1.2mm
1	HR240	Curved Strabismus Scissors 11sm
1	HR941	Skin Hook Standard 3mm
1	HR300	No3 BP Handle
1	20-7002	Scissors Iris Straight 10cm Sharp Tip
10		Swab Non Woven 10x10cm
1		Aperture Drape
1		Hand Towel
1		Tray with Gallipot Holder
1		Wrap

### 074 - Excision Set



Qty	Code	Description
1	HR121	Forceps Adson Fine Toothed
1	HR521	Needles Holder Halsey Smooth
1	HR301	Barron BP Handle
1	HR940	Skin Hook Fine 2mm
1	HR283	Scissors Tenotomy Curved Smooth Jaw
1	HR400	Forceps Mosquito Curved 13cm
10		Swab Non Woven 10x10cm
1		Aperture Drape
1		Hand Towel
1		Tray with Gallipot Holder
1		Wrap

### 039 - Excision Set Full



Qty	Code	Description
1	HR121	Forceps Adson Fine Toothed
1	HR125	Forceps Adson Non Toothed 1.5mm
1	HR300	No3 BP Handle
1	HR940	Skin Hook Fine 2mm
1	HR521	Needle Holder Halsey Smooth 13cm
1	HR240	Curved Strabismus Scissors 11sm
1	HR250	Scissors Dressing Stitch Sharp/Blunt
1	20-7002	Scissors Iris Straight Sharp 10cm
1	10-1006	Forceps Mosquito Straight
1	10-1302	Forceps Utility Lister
5		Swab Non Woven 10x10cm
1		Aperture Drape
1		60ml Gallipot
1		Hand Towel
1		Tray with Gallipot Holder
1		Wrap



# General Surgery

## Standard Sets

038 - Minor Dermatology Set



Qty	Code	Description
1	HR519	Halsey Needle Holder Smooth Narrow Jaw
1	HR300	No3 BP Handle
1	10-1302	Forceps Utility Lister
1	20-7002	Scissors Iris Straight Sharp 10cm
1	10-1006	Forceps Mosquito Straight
5		Swab Non Woven 10x10cm
1		Aperture Drape
1		60ml Gallipot
1		Hand Towel
1		Tray with Gallipot Holder
1		Wrap

020 - Minor Ops Set



Qty	Code	Description
1	HR121	Forceps Adson Fine Toothed
1	HR125	Forceps Adson Non Toothed 1.5mm
1	HR301	Barron BP Handle
1	HR523	Needle Holder Neiverts
1	20-7002	Scissors Iris Straight Sharp 10cm
5		Swab Non Woven 10x10cm
1		Drape Sheet 50x60cm
1		Hand Towel
1		Tray with Gallipot Holder
1		Wrap

047 - Minor Ops Set Full



Qty	Code	Description
1	HR122	Forceps Adson Toothed 1.2mm
1	HR940	Skin Hook Fine 2mm
1	HR283	Scissors Tenotomy Curved Smooth Jaw
1	HR400	Forceps Mosquito Curved 13cm
1	HR300	No3 BP Handle
1	HR521	Needles Holder Halsey Smooth
1	20-7002	Scissors Iris Straight Sharp 10cm
1	10-1302	Forceps Utility Lister
10		Swab Non Woven 10x10cm
1		Aperture Drape
1		60ml Gallipot
1		Tray with Gallipot Holder
1		Wrap

263 - Carpal Tunnel Set



Qty	Code	Description
1	HR122	Forceps Adson Toothed 1.2mm
1	HR301	Barron BP Handle
1	HR521	Needle Holder Halsey Smooth
1	HR240	Curved Strabismus Scissors 11mm
1	HR703	Dissector Mcdonald
1	HR951	Ramplsey Sponge Holder
1	HR283	Scissors Tenotomy Curved Smooth Jaw
2	HR700	Retractor Kilner Catpaw
5		Swab Non Woven 5x5cm
1		Tray with Gallipot Holder
1		Wrap

# General Surgery

## Forceps

**HR100 - Forceps Angled at 45° Jewellers**  
3mm Tip O/L 110mm

**HR101 - Forceps Angled at 45° Jewellers**  
7mm Tip O/L 110mm



**HR110 - Forceps Gillies Toothed**  
21mm Tying Platform O/L 110mm



**HR111 - Forceps De Bakky Toothed**  
20mm Tying Platform O/L 150mm



**HR120 - Forceps Adson Toothed**  
Std Tip Dissecting O/L 120mm


**HR122 - Forceps Adson Toothed**  
Fine Tip Dissecting O/L 120mm



**HR125 - Forceps Adson Non-Toothed**  
O/L 120mm

**HR126 - Forceps Adson Brown**  
O/L 120mm

**HR127 - Forceps Assistant Delicate**  
11mm Tying Platform O/L 115mm





# General Surgery

## Forceps


**HR130 - Forceps Iris 11cm Serrated**  
15mm Tying Platform O/L 113mm



**HR137 - Forceps Epilation**  
Angled Tip O/L 91mm



**HR160 - Forceps Allis Tissue Tying 15cm**  
5:4 Teeth O/L 150mm



**HR150 - Forceps McIndoe 15cm**  
18mm Platform O/L 155mm



**HR180 Forceps Braun Tenaculum**  
O/L 250mm



**HR183 - Forceps Vulsellum**  
O/L 250mm



**HR200 - Scissors Iris Straight**

Sharp Tip O/L 100mm



**HR210 - Scissors Iris Curved**

Sharp Tip O/L 100mm



**HR240 - Scissors Strabismus Curved**

Rounded Tip O/L 115mm



**HR241 - Scissors Strabismus Straight**

Rounded Tip O/L 115mm



**HR250 - Scissors Dressing Stitch 15cm**

Sharp/Blunt Tip

**HR251 - Scissors Dressing Stitch 13cm**

Sharp/Blunt Tip



**HR260 - Scissors Mayo Curved**

Rounded Tip O/L 135mm

**HR261 - Scissors Mayo Curved**

Rounded Tip O/L 185mm





# General Surgery

## Scissors

**HR262 - Scissors Mayo Straight**  
Rounded Tip O/L 140mm



**HR270 - Scissors Metzembaum Curved 18cm (McIndoe)**  
Rounded Tip O/L 175mm



**HR271 - Scissors Metzembaum Straight**  
Rounded Tip O/L 180mm



**HR272 - Scissors Aufricht Curved**  
Curved Sharp Tip O/L 160mm



**HR280 - Scissors Straight**  
Sharp Tip O/L 125mm



**HR281 - Scissors Stevens Tenotomy Straight**  
Sharp Tip O/L 110mm  
**HR285 - Scissors Stevens Tenotomy Straight**  
Blunt Tip O/L 115mm



**HR282 - Scissors Stevens Tenotomy Straight**  
Supercut O/L 125mm



**HR283 - Scissors Stevens Tenotomy Straight**  
Sharp Tip O/L 106mm  
**HR286 - Scissors Stevens Tenotomy Curved**  
Blunt Tip O/L 115mm



**HR300 - Handle Scalpel BP No.3**  
Rounded End O/L 125mm



**HR301 - Handle Barron BP No.3**  
O/L 135mm



**HR400 - Mosquito Curved**  
O/L 130mm

**HR420 - Mosquito Straight**  
O/L 130mm



**HR405 - Mosquito Curved**  
O/L 90mm

**HR410 - Mosquito Straight**  
O/L 90mm



**HR424 Forceps Spencer Wells**  
O/L 175mm



**HR525 Needle Holder Derf**  
O/L 125mm





# General Surgery

## Needle Holders

**HR500 - Needle Holder Mayo**  
O/L 150mm



**HR526 - Needle Holder Parkhouse**  
Angled Thumb Grip O/L 130mm



**HR510 - Needle Holder Crilewood**  
O/L 130mm



**HR519 - Needle Holder Halsey**  
Smooth / Fine O/L 130mm



**HR527 - Needle Holder Halsey**  
Smooth Angled Thumb Grip O/L 130mm



**HR520 - Needle Holder Halsey**  
O/L 130mm



**HR521 - Needle Holder Halsey**  
Smooth O/L 130mm



**HR522 - Needle Holder Kilner**  
O/L 130mm



**HR523 - Needle Holder Neivert**  
O/L 130mm



The **Blink Medical** BiPolar range has now been increased to reflect the popularity of this instrument. The forceps which come complete with their leads offer a cost effective disposable alternative to the traditional reusable instrument. All forceps in the range are produced from the highest quality components and each one goes through multiple inspections here in the UK along with a comprehensive electrical test which ensures complete safety and efficacy.



### HR630 - BiPolar McPherson

Straight with Lead O/L 110mm with universal connectors



### HR600 - BiPolar McPherson

Straight with Lead O/L 110mm



### HR601 - BiPolar McPherson

Straight with Lead O/L 200mm



### HR602 - BiPolar Adson Straight

with Lead O/L 150mm



### HR603 - BiPolar Bayonet

with Lead O/L 190mm





# General Surgery

## Retractors

**HR700 - Retractor Catspaw Kilner**  
8mm Claw Diameter O/L 150mm




**HR703 - Dissector McDonald**  
O/L 195mm



**HR704 - Retractor self retaining**  
West (Sharp) O/L 135mm

**HR708 - Retractor self retaining**  
West (Blunt) 2:3 Teeth O/L 130mm



**HR776 - Hook Cawthorne**  
O/L 150mm



**HR777 - Retractor Catspaw Mini**  
5mm Claw Diameter O/L 165mm



**HR760 Formby Wax Hook**  
O/L 170mm



**HR939 - St Bart's Hook**  
O/L 170mm



**HR820 - Towel Clip Ball & Socket**  
O/L 100mm



**HR800 - Towel Clip Backhaus**  
O/L 95mm



**HR949 - Sponge Holder Rampley**  
O/L 240mm

**HR951 - Sponge Holder Rampley**  
O/L 180mm



**HR181 - Forceps Bonney Polypus**  
O/L 25cm



**HR940 - Skin Hook Fine**  
2mm Hook O/L 165mm

**HR941 - Skin Hook Standard**  
3mm Hook O/L 165mm

**HR942 - Skin Hook McIndoe**  
5mm Hook O/L 165mm

**HR775 Skin Hook Double**  
5mm Hook

**HR944 Phlebectomy Hook**





# Cosmetic Surgery

## Forceps

### HR100 - Forcep Angled Jewellers

3mm Tip O/L 110mm

### HR101 - Forcep Angled at 45° Jewellers

7mm Tip O/L 110mm



### HR140 - Forcep Straight

Serrated Tip 8mm O/L 95mm



### HR141 - Forcep Angled Serrated Cutting

Serrated Tip 7mm O/L 95mm



### HR142 - Forcep Jewellers No.5

O/L 110mm



### HR143 - Forcep Watch Maker No.7

13mm Tip O/L 115mm



**HR710 - Fine Crocodile**  
6mm x 1.3mm Jaw



**HR711 - Forceps Tilley Nasal Long**  
18mm Serrated Jaw



**HR712 - Forceps Tilley Nasal Short**  
12mm Serrated Jaw



**HR714 - Forceps Crocodile**



**HR716 - Forceps Extra Fine Crocodile, Black Coating**  
6mm x 1.2 Jaw



**HR715 - Forceps Extra Fine Crocodile**  
5mm x 1.2 Jaw



**HR761 - Forceps Troleitsch**  
O/L 120mm



**HR722 - Ear Speculum 2mm**  
Plastic non sterile bag 100

**HR723 - Ear Speculum 3mm**  
Plastic non sterile bag 100

**HR724 - Ear Speculum 4mm**  
Plastic non sterile bag 100

**HR725 - Ear Speculum 5mm**  
Plastic non sterile bag 100

**HR726 - Ear Speculum 6mm**  
Plastic non sterile bag 100

**HR727 - Ear Speculum 7mm**  
Plastic non sterile bag 100



**T013 - Ear Speculum 2mm**

**T014 - Ear Speculum 3mm**

**T015 - Ear Speculum 4mm**

**T016 - Ear Speculum 5mm**

**T017 - Ear Speculum 6mm**

**T018 - Ear Speculum 7mm**



**HR740 - Tongue Depressor**



**HR730 - Thudichum Small**

**HR731 - Thudichum Medium**

**HR732 - Thudichum Large**



**HR751 - Jobson Horne Probe**  
125mm O/L

**HR752 - Jobson Horne Probe**  
185mm O/L



**HR763 - Formby Cerumen Scoop and Hook**



**HR830 - Curette Volkmann 5mm / 6mm cups**





**HR753 - Lempert Suction**  
2mm O/D



**HR755 - Yankauer Suction Tube**



**HR756 - Fraiser Suction Tube**  
4mm O/D



**HR750 - Zoellner Suction Tube**  
Includes 20g Tip and Trocar



**HR950 - Dental Syringe**



**HR720 - ENT Mirror**  
25mm Mirror Diameter



## HR128 - Forceps Baby Tischler

Shaft Length 210mm Jaw size 7 x 3 x 3mm



## HR129 - Forceps Eppendorfer Biopsy

Shaft length 230mm Jaw Size 7 x 3.5 x 4mm



## HR134 - Forceps Tischler Morgan Biopsy

Shaft Length 205mm Jaw Size 6 x 5 x 3mm



## HR147 - Ultimate Eppendorfer

Shaft length 235mm Jaw Size 7 x 5 x 2mm



## HR132 Forceps Vasectomy

4mm Dia Loop



## HR702 Retractor Lagenbeck



[illegible]



## Notes

[illegible]

To learn more about Corza Medical and  
how we can support you, contact us:

**[salesuk@corza.com](mailto:salesuk@corza.com) | 0121 386 8433**

Blink Medical Limited t/a Corza Medical,  
Radway Road, Solihull, West Midlands B90 4NS

