

**blink** SINGLE-USE  
INSTRUMENTS

# Single-use Ophthalmic Instruments



The Blink Medical brand was established in 2004 when Roger Tyler, the founder of Blink Medical, recognised the future needs of ophthalmic hospitals and clinics for premium, quality, single-use ophthalmic instruments.

Since this time, Blink Medical has become widely recognised as a market-leading company in the manufacture and supply of single-use ophthalmic instruments and custom packs.

Blink single-use instruments are sold globally by Corza Medical, alongside premium reusable ophthalmic instrument brands such as Katena and ASICO. We provide our customers with solutions to meet their ever-changing needs.

Corza Ophthalmology is a leading provider of single-use and reusable instruments and devices. Our portfolio includes established brands such as Blink Medical, Sharpoint, Katena, Barron and EagleVision.

Corza Medical is a growing global medical device and biosurgery company comprising innovative surgical technologies. We provide products for Ophthalmology, Biosurgery, Wound Closure and Surgical Dermatology.



Tying Forceps	4
General Forceps	4
Fixation Forceps	5
Jewellers Forceps	6
Mosquito Forceps	6
Capsulorhexis Forceps	7
Sprung Scissors	8
Iris Scissors	8
Tenotomy Scissors	8
Speculums	9
Manipulators	10
Choppers	10
Double-ended Manipulators/Choppers	11
Needle Holders	11
Refractive	12
Corneal	12
Markers	12
Sub Tenon Anaesthesia	13
Oculoplastics	14



## Manufacturing and quality

All Blink single-use instruments are manufactured to the highest standards by ISO 13485:2016 accredited manufacturing facilities. Following the initial stages of manufacturing, all products are subject to strict inspection criteria at our facility in the UK to ensure they meet the exacting standards expected by our customers. Following quality control, all items are packaged in our validated Class 8 cleanroom prior to sterilisation. All Blink single-use instruments are fully compliant with all regulatory requirements for the markets in which they are made available.



## Customer service and distribution

Our UK Customer Services and Logistics Hub in Solihull, UK, handles both Corza Medical direct distribution for the UK and international distribution for Blink single-use instruments and custom packs. Processing thousands of orders and shipments each year, our team is committed to handling your orders in a professional manner and ensuring all deliveries reach you on time.



# Forceps

## Tying Forceps

Tying forceps are available in five different designs, providing clinicians with a choice of straight and angled tips in both standard and short handle designs. The tying platforms are precisely manufactured to provide maximum grip and have smooth edges to prevent damage to sutures.

- **Tying Forceps Straight**
- 10-9009 | 5mm Tying Platform, O/L 90mm
- 10-9008 | 6mm Tying Platform, O/L 110mm
  
- **Tying Forceps Angled**
- 11-6016 | 5mm Platform Angled, O/L 90mm
- 11-6015 | 5mm Platform Angled, O/L 105mm
  
- **Tying Forceps Curved**
- 11-6012 | 5mm Platform Angled, O/L 108mm



## General Forceps


- **Max Fine Forceps Curved**
- 11-6009 | Smooth Pointed Tip, O/L 110mm
  
- **Moorfields Forceps**
- 10-1301 | 16mm Serrated Platform, O/L 107mm
  
- **Moorfields Suture Utility Forceps**
- 10-1311 | 14mm Narrow Serrated Platform, O/L 107mm
  
- **Epilation Forceps**
- 10-1300 | Angled Tip, O/L 91mm



## Fixation Forceps

Multiple handle designs are available with both notched and toothed tips. All tip designs are manufactured to provide a high level of grip for stabilising the eye during surgery whilst minimising trauma.

- **Hoskins Forceps Straight**
- 10-8008 | 0.3mm Notched, O/L 90mm
- 10-8009 | 0.3mm Notched, O/L 108mm
  
- **Pierce Colibri Forceps**
- 11-7009 | 0.3mm Notched, O/L 84mm
  
- **Bonn Forceps Straight**
- 10-1106 | 0.25mm Toothed, O/L 106mm
- 10-1107 | 0.25mm Toothed, O/L 90mm
  
- **Colibri Forceps**
- 11-7010 | 0.25mm Toothed, O/L 84mm
- 11-7011 | 0.25mm Toothed, O/L 108mm (108mm version shown)
  
- **Lister Utility Forceps**
- 10-1302 | 0.25mm Toothed, O/L 84mm
  
- **St Martins Forceps**
- 10-1203 | 0.6mm Toothed, O/L 90mm
- 10-1204 | 0.3mm Toothed, O/L 90mm





# Forceps

## Jewellers Forceps

Available in five different designs, jewellers forceps provide clinicians with a choice for holding and manipulating fine tissues that can be easily damaged by other instruments.



**Jewellers Forceps #1**  
10-2001



**Jewellers Forceps #3**  
10-2003



**Jewellers Forceps #5**  
10-2005



**Jewellers Forceps #5 Angled**  
11-2005



**Jewellers Forceps #7**  
11-2007

## Mosquito Forceps

Mosquito forceps are available in both straight and curved options for grasping and mobilising small vessels during oculoplastic surgery.



**Mosquito Forceps**  
11-1005 | Curved, O/L 100mm



**Mosquito Forceps**  
10-1006 | Straight, O/L 100mm

## Capsulorhexis Forceps




### Capsulorhexis Forceps Ultrata

- 11-5013 | Round Handle, Curved Shaft, Cystotome Tip, O/L 105mm
- 11-5012 | Round Handle, Straight Shaft, Cystotome Tip, O/L 105mm
- 11-5014 | Round Handle, Curved Shaft, Grasping Tip, O/L 105mm
- 11-5015 | Round Handle, Straight Shaft, Grasping Tip, O/L 105mm




### Capsulorhexis Forceps Ultrata

- 11-5050 | Flat Handle, Straight Shaft, Cystotome Tip, O/L 105mm
- 11-5051 | Flat Handle, Curved Shaft, Cystotome Tip, O/L 105mm
- 11-5052 | Flat Handle, Straight Shaft, Grasping Tip, O/L 105mm
- 11-5053 | Flat Handle, Curved Shaft, Grasping Tip, O/L 105mm



#### Cystotome Tip

A sharp tip design that can be used for creating the initial flap and grasping the lens capsule when carrying out a capsulorhexis.



#### Grasping Tip

Designed with larger grasping platforms for carrying out a capsulorhexis when the initial flap has been created with another instrument.

### Capsulorhexis Forceps Inamura

- 11-5021 | Flat Handle, Curved Shaft, Cystotome Tip, O/L 105mm
- 11-5011 | Flat Handle, Straight Shaft, Cystotome Tip, O/L 105mm



### Capsulorhexis Forceps Inamura

- 11-5020 | Round Handle, Curved Shaft, Cystotome Tip, O/L 105mm


# Scissors | Speculums

## Sprung Scissors

Sprung scissors are available in Westcott, Vannas and Ong designs. The Westcott designs have a variety of tip forms, as well as being available in straight or curved blade options. The drilled handle design increases grip and reduces the weight of the scissors, resulting in improved tactility.

### Westcott Scissors

- 20-4013 | Straight, O/L 120mm
- 20-4015 | Straight, Sharp Tip, 120mm
- 21-4010 | Curved, O/L 120mm
- 21-4014 | Curved, Sharp Tip, O/L 120mm
- 21-4011 | Curved, Dolphin Nose, O/L 120mm
- 21-5001 | Curved, O/L 100mm



### Vannas Scissors

- 20-6000 | Straight, O/L 90mm
- 21-6001 | Curved, O/L 90mm

### Ong Scissors

- 21-3010 | Angled Tip 9mm, O/L 92mm




## Iris Scissors and Tenotomy Scissors

- 20-7001 | Iris, Straight, Blunt Tip, O/L 100mm
- 20-7002 | Iris, Straight, Sharp Tip, O/L 100mm
- 21-7004 | Iris, Curved, Sharp Tip, O/L 100mm
- 21-7005 | Tenotomy, Curved, Pointed Tip, O/L 115mm
- 21-7006 | Tenotomy, Curved, Blunt Tip, O/L 115mm


## Speculums

Speculums provide access to the eye for all types of ophthalmic procedures. The Blink speculum range includes a wide variety of designs. Most popular are Barraquer wire speculums, which are available in eight different designs and offer multiple benefits during ophthalmic surgery. Open Loop are widely used for cataract surgery, whereas Closed Loop and Solid Blade are favoured for administering ocular injections.



### Barraquer Closed Loop

- 30-2004 | Straight, Width 30mm, Length 43mm
- 31-2003 | Temporal, Width 30mm, Length 42mm
- 31-2002 | Nasal, Width 30mm, Length 42mm



### Kratz-Barraquer Open Loop

- 31-1006 | Temporal, Width 30mm, Length 40mm
- 30-1007 | Straight, Width 30mm, Length 42mm
- 31-1005 | Nasal, Width 30mm, Length 40mm
- 31-1009 | Temporal, Curved, Width 30mm, Length 42mm



### Leibinger Speculum

- 31-3002 | Wire Loop, Temporal, O/L 65mm



### Leibinger Speculum

- 31-3000 | Solid Blade, Temporal, O/L 65mm



### Barraquer Speculum Solid Blade

- 30-4000 | Temporal, Width 30mm, Length 44mm



### Clarkes Speculum

- 31-3004



### Burratto Speculum

- 30-3001 | Solid Blade, O/L40mm



### Cook Speculum

- 31-3003 | Solid Blade, O/L40mm



### Paediatric Speculum Spring

- 30-3002



### Neonatel Speculum Angled Temporal

- 31-3100



### Paediatric Speculum

- 31-3102

# Manipulators | Needle Holders



## Manipulators

Available in ten different designs, the Blink manipulator range provides clinicians with a wide choice of options.



**Sinskey Lens Hook Straight**  
40-5002



**Sinskey Lens Hook Angled**  
41-5001



**Nucleus Y Shaped Angled Horizontal Axis**  
41-2005



**Nucleus Y Shaped Angled Vertical Axis**  
41-2006



**Spatula Manipulator**  
41-3004



**Kuglen Hook Angled**  
41-4014



**Mushroom Tipped Manipulator**  
41-8009



**IOL Lester Manipulator**  
41-1007



**Spatula Neuhann Manipulator**  
41-3005



**Connor Wand**  
41-4013

## Choppers

Available in eight different designs to meet clinicians' requirements when carrying out lens removal, these designs provide options for vertical and horizontal chopping techniques.



**Drysdale Paddle Chopper**  
41-6007



**Sadiq Chopper**  
41-4005



**Koch Chopper 1.4mm**  
41-4008



**Rosen Chopper 1.2mm**  
41-4010



**Seibel Chopper**  
41-4012



**Nagahara Chopper**  
41-4009



**Greenes Chopper**  
41-7008



**McPhearson Chopper**  
41-3006

## Double-ended Manipulators/Choppers

By bringing together our most popular designs in this range, these four options of double-ended manipulators/choppers provide the benefits of two products in one.



**Sadiq Chopper/Mushroom Tipped Manipulator**  
41-9012



**Sinskey Manipulator/Mushroom Tipped Manipulator**  
41-9010



**Rosen Chopper/Mushroom Tipped Manipulator**  
41-9013



**Nagahara Chopper/Connor Wand**  
41-9014

## Needle Holders



**Needle Holder Castroviejo**  
90-7001 | 12mm Tying Platform, Straight With Lock, O/L 135mm

**Needle Holder Castroviejo**  
90-7002 | 12mm Tying Platform, Straight, O/L 135mm

**Needle Holder Castroviejo**  
91-7003 | 12mm Tying Platform, Delicate Curved, O/L 135mm

**Needle Holder Troutman**  
91-7004 | 6mm Tying Platform, Curved, O/L 105mm



## Refractive



• **Toric Alignment Marker, 3, 6, 9 O'Clock**  
• 91-6001



• **Toric Degree Gauge 11.7mm I.D**  
• 91-6000



• **Toric IOL Marker 11mm Wide**  
• 91-6002

## Corneal

• **Epithelial Spatula**  
• 51-0001

• **Epithelial Spatula, Rounded Tip**  
• 51-0003

• **Alcohol Well**  
• 51-0002 | 9.5mm



## Markers

• **Double-ended Marker**  
90-1015 | Flat Handle, 3.5mm-4mm, O/L 113mm

• **Braunstein Marker**  
90-1017 | Round Handle, 3.5mm-4mm, O/L 80mm

• **Intravitreal Caliper/Fixator**  
90-1016 | Round Handle, 3.5mm, O/L 70mm

• **Caliper & Rule**  
90-1014 | 0-20mm Measurement Span

## Sub Tenon Anaesthesia

Corza Medical is focussed on providing healthcare professionals with innovative new products that save time, money and effort.

One such product is our Conjunctival Probe which is uniquely designed to create the perfect-sized incision for safe and effective injection of anaesthetic during the Sub Tenon's procedure. The tip has been designed to be the same size as the Sub Tenon Cannula, resulting in far less backflow of anaesthetic out of the wound, reduced red eye due to accidental puncturing of blood vessels, and improved akinesia due to the more effective distribution of anaesthetic. The SISTA technique (Small Incision Sub Tenon Anaesthesia) can be viewed below.

Should one of our standard sets not meet your requirements, please contact our helpful customer service staff who would be happy to help you specify your own custom pack.



• **Conjunctival Probe**  
• CP01



Raise the conjunctiva using a pair of non-toothed Moorfields forceps, holding them at 90 degrees to the direction of the probe.




Holding the probe parallel to the surface of the sclera, puncture the conjunctiva.



Keeping the probe in contact with the sclera, follow the curve of the globe for approximately 1cm.



Withdraw the probe but maintain the grip on the conjunctiva with the forceps.



Insert the Sub Tenon Cannula to inject the anaesthetic in the usual manner.



• **001 Sub Tenon Instrument Set**  
• Scissors Westcott Curved Standard Tip  
• Speculum Kratz-Barraquer Angled Temporal  
• Moorfields Forceps



• **270 CP Sub Tenon Set**  
• Conjunctival Probe  
• Speculum Kratz-Barraquer Angled Temporal  
• Moorfields Forceps



• **300 Complete Sub Tenon Set**  
• Scissors Westcott Curved Fine Tip  
• Speculum Barraquer Angled Temporal  
• Moorfields Forceps  
• Sub Tenon Cannula  
• Syringe 5ml Luer Lock  
• Swab Non Woven 10x10cm  
• Tray with Gallipot Holder



• **290 Complete CP Sub Tenon Set**  
• Speculum Barraquer Angled temporal  
• Moorfields Forceps  
• Conjunctival Probe  
• Syringe 5ml Luer Lock  
• Swab Non Woven 10x10cm  
• Tray with Gallipot Holder

# Oculoplastics



- **Liebreich Probes**
- **90-2007** | Liebreich Probe 0000+000 | O/L 125mm
- **90-2008** | Liebreich Probe 000+00 | O/L 125mm
- **90-2009** | Liebreich Probe 00+0 | O/L 125mm
- **90-2010** | Liebreich Probe 1+2 | O/L 125mm
- **90-2011** | Liebreich Probe 3+4 | O/L 125mm

## Curettes

- **Curette Double-ended**
- **90-3007** | 2mm & 3mm Cups, O/L 130mm
- **Curette Double-ended**
- **90-3008** | 2mm & 5mm Cups, O/L 130mm
- **Curette Double-ended**
- **90-3009** | 2mm & 4mm Cups, O/L 130mm



- **Clamp Chalazion**
- **90-4005** | Small, 11mm Diameter Clamp, O/L 90mm
- **90-4006** | Medium, 18mm x 11mm Clamp, O/L 80mm
- **90-4007** | Large, 23mm x 15mm Clamp, O/L 90mm



- **Desmares Entropion Clamp**
- **90-4008** | Left, Clamp Diameter 24mm x 16mm, O/L 95mm



- **Desmares Entropion Clamp**
- **90-4009** | Right, Clamp Diameter 24mm x 16mm, O/L 95mm



- **Desmarres Retractor**
- **90-9013** | 14mm Wide Blade, O/L 130mm



- **Graefe Strabismus Hook**
- **90-5004** | 10mm Tip Length, O/L 140mm



- **Alexandra Muscle Hook**
- **90-5005** | 11mm Tip Length, O/L 130mm



- **Dilator Punctal Double-ended**
- **90-1011** | Pointed Tips, O/L 138mm



- **Dilator Nettle Ship**
- **90-1012** | 20mm Tapered Tip, O/L 118mm




- **Bulldog Clamp**
- **90-8012** | O/L 52mm



- **092 – Chalazion Set – Small Clamp**
- **90-4005** Chalazion Clamp Small
- **90-3009** Curette Double Ended Small
- **No.11 Scalpel**
- **Wrap 25x25cm**
- **III – Chalazion Set – Medium Clamp**
- **90-4005** Chalazion Clamp Small
- **90-3009** Curette Double Ended Medium
- **No.11 Scalpel**
- **Wrap 25x25cm**

## General

- **Thornton Ring**
- **90-1020** | 13mm Diameter





To learn more about Corza Medical and how we can support you, contact us:

**[salesuk@corza.com](mailto:salesuk@corza.com) | 0121 386 8433**

Blink Medical Limited t/a Corza Medical,  
Radway Road, Solihull, West Midlands B90 4NS  
COR-ADV-12122022-Rev1

 **corzamedical**  
Ophthalmology