

**EagleVision<sup>®</sup>** PUNCTUM  
PLUGS

# Punctum Plugs



EagleVision was founded by Jerre Freeman, MD, in 1983. He invented and patented the world's first Punctum Plug nine years earlier, in 1974. Driven by excellence, EagleVision has become the brand that is synonymous with punctum plugs across the globe.

The first tapered shaft plugs were developed in 1994 and the first low profile rim plugs in 1998. Since this time EagleVision has continued to invest in new designs to provide improved patient comfort and fit, such as SuperEagle and SuperFlex.

EagleVision was acquired by Katena Products in 2016, bringing together two of the most trusted brand names in the ophthalmic industry. Now, under Corza Medical, EagleVision continues to be the most trusted brand name in the global punctum plug market.

Corza Ophthalmology is a leading provider of single-use and reusable instruments and devices. Our portfolio includes established brands such as Blink Medical, Sharpoint, Katena, Barron and EagleVision.

Corza Medical is a growing global medical device and biosurgery company comprising innovative surgical technologies, we provide products for Ophthalmology, Biosurgery, Wound Closure and Surgical Dermatology.



Available in boxes of 2

## Plug 1

Unique dual-lobed design allows the plug to fit a wide range of Punctum sizes

Manufactured from a soft grade of medical silicone to enhance patient comfort

Ideal option for users with very low and intermittent demand for Punctum Plugs

Individually sterile packaged, preloaded on a disposable inserter/dilator for ease of use

## SuperEagle

Tapered shaft and pivoting wide-flex nose allow for better retention

Individually sterile packaged, preloaded on a disposable inserter/dilator for ease of use

Available in three sizes for convenience

- 0.4mm – 0.6mm
- 0.6mm – 0.8mm
- 0.8mm and larger

Available in boxes of 2



Available in boxes of 2

## SuperFlex

Contouring traction ribs provide ultimate flexibility for easy insertion and superior retention

Low-profile rim and soft, low durometer silicone enhances patient comfort

Available in ten different sizes, from 0.3mm to 1.2mm diameter, to provide the perfect fit

Individually sterile packaged, preloaded on a disposable inserter/dilator for ease of use



Corza Code	Description	Size	Unit of Purchase
30012	Plug 1	One size to fit all	Box of 2
31302	Punctum Plug - SuperEagle	Small: 0.4mm – 0.6mm	Box of 2
31312	Punctum Plug - SuperEagle	Medium: 0.6mm – 0.8mm	Box of 2
31322	Punctum Plug - SuperEagle	Large: 0.8mm and larger	Box of 2
31732	Punctum Plug - SuperFlex	0.3mm x 1.1mm	Box of 2
31742	Punctum Plug - SuperFlex	0.4mm x 1.1mm	Box of 2
31752	Punctum Plug - SuperFlex	0.5mm x 1.2mm	Box of 2
31762	Punctum Plug - SuperFlex	0.6mm x 1.3mm	Box of 2
31772	Punctum Plug - SuperFlex	0.7mm x 1.4mm	Box of 2
31782	Punctum Plug - SuperFlex	0.8mm x 1.5mm	Box of 2
31792	Punctum Plug - SuperFlex	0.9mm x 1.7mm	Box of 2
31702	Punctum Plug - SuperFlex	1.0mm x 1.8mm	Box of 2
31712	Punctum Plug - SuperFlex	1.1mm x 2.0mm	Box of 2
31722	Punctum Plug - SuperFlex	1.2mm x 2.0mm	Box of 2
31802	Punctum Plug - SuperFlex	1.3mm x 2.0mm	Box of 2
2040	Coroneo Punctal Gauge - Reusable		Each

To learn more about Corza Medical and how we can support you, contact us:

**[salesuk@corza.com](mailto:salesuk@corza.com) | 0121 386 8433**

Distributed by:  
Blink Medical Limited t/a Corza Medical  
Radway Road, Solihull, West Midlands, B90 4NS, UK

Manufactured by:  
Katena Products inc (a Corza Medical Company)  
6 Campus Drive, Suite 310, Parsipanny, New Jersey, 07054, USA

BLNK-CAT1-11282022-Rev0

 **corzamedical**  
Ophthalmology