

**SharpPoint** MICROSURGICAL  
KNIVES

**Single-use  
Ophthalmic Knives**



# A rich history empowering a bright future.

The Sharpoint story began in 1969 when brothers William and Rolf Schmidt co-founded Sharpoint Inc. From their garage workshop, they developed the technology that would set the standard for producing ultra-sharp, premium ophthalmic surgical knife and blade products for the next 50 years. Now, under the ownership of Corza Medical, Sharpoint knives and blades are globally recognised by Ophthalmic surgeons for their sharpness, consistency and varying designs that meet the requirements of surgeons today. With continued investment in new products such as the CC3D+ Slit Knives, Corza Medical is committed to providing surgeons globally with the best choices in knives and blades for the next 50 years.

Corza Ophthalmology is a leading provider of single-use and reusable instruments and devices from well-known brands such as Blink Medical, Sharpoint, Katena, Barron and EagleVision. Corza Ophthalmology is proud to offer a range of proven solutions designed and manufactured to help surgeons be at their best at every step, in every procedure. Corza Medical is a growing global medical device and biosurgery company of innovative surgical technologies delivering products for Ophthalmology, Biosurgery, Wound Closure and Surgical Dermatology. Our mission is to champion and aid surgical providers with remarkable service and exceptional choices that save time, money and effort. We promise to always prioritise you, so that you can prioritise patients.




## Contents

<b>Sharpoint Safety Knives</b>	
Safety Knives .....	01
<b>Phaco/Implant Incision Knives</b>	
Clear Control 3D Plus Slit Knives .....	02
Slit Knives, Double-Bevel .....	02
Slit Knives, Bevel-Up .....	02
ClearTrap™ Trapezoid Clear Corneal/Implant Knives, Bevel-Up .....	03
Pilot Tip™ Implant Knives, Bevel-Up .....	03
<b>Paracentesis Incision Knives</b>	
ClearPort™ Paracentesis Knives .....	04
Stab Knives .....	04
<b>Vitrectomy Knives</b>	
Vitrectomy Knives .....	05
<b>Crescent Knives</b>	
Sharptome™ Crescent Knives .....	06
Mini-Crescent Knives .....	06
<b>Groove and Limbal Relaxation Incision Knives</b>	
Precision/Restricted Depth Knives .....	07
<b>Specialty Knives</b>	
Spoon Knives, Double-Bevel .....	08
Specialty Knives .....	08
Blade and Tip Handle .....	08



## Pioneer in Blade Precision.

Over 50 years of blade edge design experience – Sharpoint Safety Knives, with Infinite™ Edge Technology, are designed to provide:

- Accurate incisions
- Exceptional wound architecture
- Promote self-sealing wounds



The **Infinite Edge™ Technology** of blade etching and electro surface polishing produces knives of exceptional sharpness, outstanding precision and consistency.



### ● Crescent

Designed to create a smooth scleral shelf to promote a self-sealing wound, the blade's consistent sharpness virtually eliminates tissue drag and distortion.



### ● Sideport MVR

Designed to create consistent, accurate and self-sealing paracentesis incisions.



### ● Sideport Stab

Provides a strong and true point for easy penetration during initial stab incisions.



### ● Slit Advanced Double Bevel (CC3D+)

Designed to facilitate instrument insertion through a multidimensional incision. Micro bevels add vertical dimension, placing less stress on wound margins.



### ● Slit (Single and Double-Bevel)

Sharpoint angled precision slit knives penetrate easily while creating an accurate width incision.



# SharpPoint Safety Knives

- Single-Handed "No-Look" Blade Activator
- Retractable Blade prevents sharps injuries during handling and use
- Slim No-Hood Tip design allows superior visualisation during use
- Snap Lock Slider provides tactile feedback when in locked position
- Colour Coded Blade Type for quick, accurate identification
- Clear Housing Design protects blade against inadvertent damage
- 5 per box

## Safety Slit Single Bevel



Knife	Description
62-2260G	2.2 mm, Slit Angled, Single Bevel
62-2460G	2.4 mm, Slit Angled, Single Bevel
62-2560G	2.5 mm, Slit Angled, Single Bevel
62-2760G	2.75 mm, Slit Angled, Single Bevel
62-2860G	2.8 mm, Slit Angled, Single Bevel
62-3060G	3.0 mm, Slit Angled, Single Bevel

## Safety Slit Double Bevel



Knife	Description
62-2442G	2.4 mm, Slit Angled, Double Bevel
62-2542G	2.5 mm, Slit Angled, Double Bevel
62-2742G	2.75 mm, Slit Angled, Double Bevel
62-3042G	3.0 mm, Slit Angled, Double Bevel

## Safety Crescent Knives




Knife	Description
64-1030G	2.3 mm, Crescent, 55 Degree Bevel Up

## Safety Slit Advanced Double Bevel




Knife	Description
62-2250G	2.2 mm, CC3D+ Slit, Angled, Advanced Double Bevel
62-2450G	2.4 mm, CC3D+ Slit, Angled, Advanced Double Bevel

## Safety Sideport Stab



Knife	Description
62-1501G	15 Degree, Sideport, Straight, Stab

## Safety Sideport MVR Knives



Knife	Description
68-2010G	1.0 mm, Sideport MVR, Angled


# Phaco Incision Knives

## Clear Control 3D Plus Slit Knives

SharpPoint™ Clear Control 3D Plus slit knives are designed to facilitate instrument insertion through a multidimensional incision. Micro bevels are designed to add vertical dimension, placing less stress on wound margins, allowing for instrument insertion and manipulation. 6 per box.



Knife	Description
72-1650	1.6mm, Angled, Full Handle
72-1850	1.8mm, Angled, Full Handle
72-2050	2.0mm, Angled, Full Handle
72-2250	2.2mm, Angled, Full Handle
72-2450	2.4mm, Angled, Full Handle
72-2550	2.5mm, Angled, Full Handle
72-2650	2.65mm, Angled, Full Handle
72-2750	2.75mm, Angled, Full Handle
72-2850	2.85mm, Angled, Full Handle
72-3050	3.0mm, Angled, Full Handle
72-3250	3.2mm, Angled, Full Handle



## Slit Knives, Double-Bevel

SharpPoint™ double-bevel slit knives easily create accurate width incisions with virtually no tissue drag or distortion. 6 per box.



Knife	Description
72-1831	1.8mm, Angled, Double-Bevel
72-2231	2.2mm, Angled, Double-Bevel
72-2431	2.4mm, Angled, Double-Bevel
72-2531	2.5mm, Angled, Double-Bevel
72-2631	2.65mm, Angled, Double-Bevel
72-2831	2.8mm, Angled, Double-Bevel
72-2931	3.0mm, Angled, Double-Bevel
72-3231	3.2mm, Angled, Double-Bevel

## Slit Knives, Bevel-Up

SharpPoint™ angled, bevel-up precision slit knives penetrate easily while creating an accurate width incision for optimal fit to the phacoemulsification tip. 6 per box.



Knife	Description
72-2261	2.2mm, Angled, Bevel-Up
72-2461	2.4mm, Angled, Bevel-Up
72-2561	2.5mm, Angled, Bevel-Up
72-2661	2.65mm, Angled, Bevel-Up
72-2761	2.75mm, Angled, Bevel-Up
72-2861	2.85mm, Angled, Bevel-Up
72-2961	3.0mm, Angled, Bevel-Up
72-3261	3.2mm, Angled, Bevel-Up



## Phaco/Implant Incision Knives

### Trapezoid Design

Yields phaco and implant incisions in one knife

### Consistency


Allows surgeon to enlarge the phaco incision to a consistent, precise width for IOL insertion

### Wound Integrity

Designed to maintain wound integrity by eliminating the implant knife

### Indicator Lines

Patented width indicator lines at the phaco and implant mark provide a visual reference



### ClearTrap™ Trapezoid Clear Corneal/Implant Knives, Bevel-Up

Sharpoint™ ClearTrap™ knives feature a trapezoid shape design containing patented incision width indicators visible on each blade. The unique thin blade design and consistent sharpness create incisions with virtually no tissue drag or distortion. 6 per box.



Knife	Description
75-1214	1.2 - 1.4mm, Angled, Width Indicator
75-1617	1.5 - 1.7mm, Angled, Width Indicator
75-1620	1.5 - 2.0mm, Angled, Width Indicator
75-1720	1.7 - 2.0mm, Angled, Width Indicator
75-2022	2.0 - 2.2mm, Angled, Width Indicator
75-2522	2.5 - 2.9mm, Angled, Width Indicator
75-2723	2.75 - 3.0mm, Angled, Width Indicator
75-2724	2.75 - 3.2mm, Angled, Width Indicator
75-2725	2.75 - 3.5mm, Angled, Width Indicator
75-2824	2.85 - 3.2mm, Angled, Width Indicator
75-2826	2.85 - 3.75mm, Angled, Width Indicator

### Pilot Tip™ Implant Knives, Bevel-Up

Sharpoint™ Pilot Tip™ implant knives are designed to easily follow the slit blade incision while the cutting edges create a precise opening for lens insertion. 6 per box.



Knife	Description
76-5361	3.2mm, Angled, Bevel-Up
76-5161	5.2mm, Angled, Bevel-Up
76-5561	5.5mm, Angled, Bevel-Up

## Paracentesis Incision Knives

### ClearPort™ Paracentesis Knives


Technology = Superior Wound Architecture

### Exceptional Wound Architecture

Created for consistent and accurate wound architecture to achieve self-sealing wounds

### Blade Geometry

and choice of sizes allows for precise instrument fit



### Paracentesis Knives

Sharpoint™ ClearPort™ knives are designed to create consistent, accurate and self-sealing paracentesis incisions. 6 per box.



Knife	Description
78-2310	0.8mm, Angled, Full Handle
78-2020	1.0mm, Angled, Full Handle
78-2010	1.1mm, Angled, Full Handle
78-2019	1.4mm, Angled, 1.0mm and 1.2mm Indicator Lines

### Stab Knives

Sharpoint™ stab knives are heat-tempered to provide a strong, true point for easy penetration during initial stab incisions. 6 per box.



Knife	Blade/Tip	Description
72-1501	72-1502	15° Straight, Full Handle
72-1591	-	15° Right-angled
72-2201	72-2202	22.5° Straight, Full Handle
72-3001	72-3002	30° Straight, Full Handle
72-4501	72-4502	45° Straight, Full Handle



Knife	Blade/Tip	Description
72-1531	72-1532	15° 3.0mm, Straight, Restricted Blade Depth
72-1551	72-1552	15° 5.0mm, Straight, Restricted Blade Depth
72-3036	-	30° 3.5mm, Straight, Restricted Blade Depth

**Vitreotomy Knives**

SharpPoint™ vitrectomy knives create small, precise, access incisions which virtually eliminate wound leaks around probes and other posterior segment instruments. 6 per box.



Knife	Description
71-1901	19 Gauge, Straight
71-2001	20 Gauge, Straight
71-2301	23 Gauge, Straight
71-2501	25 Gauge, Straight



Knife	Description
71-2061	20 Gauge, Angled

**Sharptome™ Crescent Knives**

SharpPoint™ Sharptome™ crescent knives are designed to create a smooth scleral shelf to promote a self-sealing wound. The blade's consistent sharpness virtually eliminates tissue drag and distortion. 6 per box.



Knife	Blade/Tip	Description
74-1000	74-5000	2.0mm, Straight
----	78-6900	Micro-Blade, Straight, Double-Bevel



Knife	Blade/Tip	Description
74-1010	74-5062	2.0mm, Angled, Bevel-Up
74-1031	-	2.0mm, Angled, Double-Bevel

**Mini-Crescent Knives**

SharpPoint™ mini-crescent knives are designed to dissect a smaller scleral shelf. The blade's consistent sharpness virtually eliminates tissue drag and distortion. 6 per box.



Knife	Description
74-6910	1.25mm, Angled, Double-Bevel

### Precision Depth Knives

SharpPoint™ precision depth knives are designed to create groove and limbal relaxing incisions. Single formed design ensures precision accuracy. Dual cutting edges allow for incisions in left or right directions. 6 per box.



Knife	Description
72-0303	300 µm Straight, Full Handle
72-3503	350 µm Straight, Full Handle
72-5503	550 µm Straight, Full Handle
72-6003	600 µm Straight, Full Handle

### Restricted Depth Knives

SharpPoint™ restricted depth knives are designed to create groove and limbal relaxing incisions. Single formed design ensures precision accuracy. Dual cutting edges allow for incisions in left or right directions. 6 per box.



Knife	Description
72-0301	300 µm Straight, Full Handle
72-0401	450 µm Straight, Full Handle

### Spoon Knives, Double-Bevel

SharpPoint™ spoon knives feature a unique 320° cutting radius and an ultra-thin blade to enhance precision cutting capabilities where multidirectional techniques are performed. 6 per box.



Knife	Description
74-2261	2.2mm, Angled, Double-Bevel
74-3061	3.0mm, Angled, Double-Bevel

### Specialty Knives

SharpPoint™ specialty knives are available in a variety of shapes and sizes for various techniques and procedures. 6 per box.



Knife	Description
78-6810	Needle Blade, Straight, 4.0mm Cutting Edge



Knife	Blade/Tip	Description
-	78-5700	Scleral Blade, Straight



Knife	Blade/Tip	Description
-	78-6600	2.3mm, Lamellar Blade, Angled, Bevel-Up



Knife	Blade/Tip	Description
-	78-6900	Micro-Blade, Straight, Double-Bevel, Sharp All Around

### Blade and Tip Handle

SharpPoint™ aluminum blade and tip handle is compatible with the threaded tips of SharpPoint™ blade and tip items. The inner threads are standardised to fit most competitive collets and blades currently available. Sold individually.



Handle	Description
92-9003	16cm Length

For more information, please email  
[sales.uk@corza.com](mailto:sales.uk@corza.com) or call 0121 386 8433

Corza Medical (Blink Medical Ltd)  
Radway Road, Solihull, West Midlands, B90 4NS  
Phone: 0121 386 8433 Email: [sales.uk@corza.com](mailto:sales.uk@corza.com) [www.blinkmedical.co.uk](http://www.blinkmedical.co.uk)  
All trademarks are property of their respective owners. ©2022 Surgical Specialties Corporation  
All Rights Reserved. COR-BRO-04142022-Rev1

 **corzamedical**  
Ophthalmology