



Needle & Trocar Compatibility Guide

A guide for surgical needles that may be compatible with commercially available trocar devices.

Taper Point Needle Name	Trocar Size			
	8mm	11mm	12mm	15mm
RB-1 1/2 Circle 17mm	●	●	●	●
SH-1 1/2 Circle 22mm	●	●	●	●
SH 1/2 Circle 26mm	●	●	●	●
CT-2 1/2 Circle 26mm	●	●	●	●
UR-6 5/8 Circle 26mm		●	●	●
CT-1 1/2 Circle 36mm			●	●
CT 1/2 Circle 40mm				●

Quill® is available through distribution or direct. For questions or a demonstration, please contact your local sales representative or contact us. For procedural videos, visit youtube.com/TheQuillDevice

Corza.com • 1.877.991.1110 • Service@corza.com

Surgical Specialties Corporation, a Corza Medical company.
247 Station Drive, Suite NE1, Westwood, MA 02090 U.S.A. All trademarks are property of their respective owners. ©2021 Corza Medical. All Rights Reserved
COR-226 R0

