

# microMend<sup>®</sup> PRO

Skin Closure Device

Distributed by  **corzamedical**

## CLINICAL PRODUCT PORTFOLIO

### Product Description

microMend<sup>®</sup> PRO is a unique new alternative to sutures, surgical glues and common adhesive strips for closing the skin surface as well as protecting the incision site afterwards.

### Regulatory Status

microMend<sup>®</sup> PRO is an FDA Class 1 medical device listed with the FDA and available for sale in the United States.

## ORDERING INFORMATION

Product Code	Description	Size	Devices Per Box	Packaging
MM12XS2	microMend PRO Small	7mm x 28mm	24	2 Devices / Pack 12 Packs / Box
MM12SM2	microMend PRO Medium	12mm x 28mm	24	2 Devices / Pack 12 Packs / Box
MM12SM4	microMend PRO Medium	12mm x 28mm	48	4 Devices / Pack 12 Packs / Box
MM12LG2	microMend PRO Large	12mm x 43mm	24	2 Devices / Pack 12 Packs / Box
MM12SU2	microMend PRO Umbilical	12mm x 26mm	24	2 Devices / Pack 12 Packs / Box
MM12W2A	microMend PRO Wide With Aperture	38mm x 28mm	24	2 Devices / Pack 12 Packs / Box
MM12W2S	microMend PRO Wide Solid Seal	38mm x 28mm	24	2 Devices / Pack 12 Packs / Box

**KitoTech**  
MEDICAL

© 2022 KitoTech Medical, Inc. All Rights Reserved. microMend<sup>®</sup> is a registered trademark of KitoTech Medical, Inc.

**Manufactured For:**  
KitoTech Medical, Inc.  
500 Yale Ave N | Seattle, WA 98109  
www.micromendskinclosure.com

**Distributed by:**  
Surgical Specialties Corporation  
(a Corza Medical Company)  
247 Station Drive, Suite NE1  
Westwood, MA 02090  
www.corza.com

### Customer Service


Phone: **(877) 991-1110**  
Email: [service@corza.com](mailto:service@corza.com)



## microMend PRO Small

**Ideal for small areas, such as the fingers, toes, ears, and nose.**


One microMend PRO Small device closes a wound approximately 0.75 cm in length. Two PRO Small devices (contained in one pouch) close a wound that is approximately 1.5 cm in length.



## microMend PRO Medium

**Broadest application for closing most wounds. It is used for closing most lacerations (face, neck, and hands).**


One microMend PRO Medium device closes a wound approximately 1.5 cm in length. Two PRO Medium devices (contained in one pouch) close a wound that is approximately 3 cm in length.



## microMend PRO Large

**Longer device with extra gripping strength for closing larger wounds on the extremities, back, chest, and abdomen. Great for closing hard to repair skin tears.**


One microMend PRO Large device closes a wound approximately 1.5 cm in length. Two PRO Large devices (contained in one pouch) close a wound that is approximately 3 cm in length.




## microMend PRO Umbilical

**Designed to close umbilical port sites associated with laparoscopic and robotic surgeries.**

One microMend PRO Umbilical device is needed to close most Umbilical Port Sites (wounds of less than 1.5 cm in length). When the Port Site is longer than 1.5 cm, two devices are usually required (contained in one pouch).



**WIDE-A**



**WIDE-S**

## microMend PRO Wide

**Three times wider for easier and faster closure of longer incisions with the same holding strength as standard microMend devices.**

Ultra-flexible and elastic making it ideal for use in joint surgeries, Cesarean sections, and other procedures in areas of increased movement and varying skin architecture.

One microMend PRO Wide closes a wound that is approximately 4 cm in length. Two (2) microMend PRO Wide devices (contained in one pouch) close a wound that is approximately 8 cm in length.

Available **with apertures (WIDE-A)** to allow drainage, or **without apertures (WIDE-S)** for a solid seal.